
ThermoFisher
S C I E N T I F I C

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

- + Admin
- OMNIC
 - Access Control
 - Ability to run OMNIC
 - File
 - Open
 - Save
 - Save As
 - Save Group
 - Mail Spectra
 - Sign File
 - Sign Spectra
 - Verify File
 - Open Configuration
 - Save Configuration As
 - Open Log
 - Print
 - Printer Setup
 - Delete Files
 - Edit
 - Collect
 - View
 - Process
 - Analyze
 - Raman
 - Report
 - Window
 - Help
 - Experiment Setup tabs
 - Library Setup tabs
 - Library Manager functions
 - Custom toolbar and menu items
 - System Policies
 - Signature Meanings

Access Control

Ability to run OMNIC

Access control settings grant or deny individual users or groups access to protected functions in the client application.

List names from: AMER

Access rights:

- Bob.Manager
- Joe.Labtech
- Lee.Scientist
- Users

< - Add

Remove ->

Names:

- Administrator
- Authenticated Users
- Bob.Manager
- Everyone
- Guest
- HelpAssistant
- Joe.Labtech
- Lee.Scientist
- None
- SUPPORT_388945a0

Select from the list of names or enter a name here:

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The screenshot displays the Security Administration application window. The left pane shows a tree view with 'System Policies' selected. The main pane shows the configuration for 'System Policies', including a description, policy groups, and a list of policy members. A list of names is visible on the right side of the main pane. A detailed list of system policies is shown on the right side of the image, with the policy 'When user logs on, run only OMNIC' highlighted in blue. Two black arrows point from the 'System Policies' folder in the left pane to the detailed list on the right, and from the 'Policy group members' list in the main pane to the detailed list on the right.

System Policies




- Prevent overwriting of files
- Include original data format when saving
- Require signature when saving spectrum
- Require signature when saving experiment
- Require signature when saving configuration
- Require signature when adding report to notebook
- Prevent signature meaning entry when signing file
- Require reason entry when changing spectrum
- Deny access to menu items not in security database
- Prevent printing, copying or cutting of unsaved data
- When user logs on, run only OMNIC**
- Automatically start screen saver after (minutes)
- Prevent changing directories when opening files
- Prevent changing directories when saving files
- Disable delete, rename and right-click in file dialog boxes
- Default configuration
- Directory for spectral data
- Directory for experiments
- Directory for log files
- Directory for quant methods
- Directory for macros
- Directory for user libraries
- Directory for report templates and notebooks
- Directory for saving spectra automatically

System Policies

- Prevent overwriting of files
- Include original data format when saving
- Require signature when saving spectrum
- Require signature when saving experiment
- Require signature when saving configuration
- Require signature when adding report to notebook
- Prevent signature meaning entry when signing file
- Require reason entry when changing spectrum
- Deny access to menu items not in security database
- Prevent printing, copying or cutting of unsaved data
- When user logs on, run only OMNIC
- Automatically start screen saver after (minutes)
- Prevent changing directories when opening files
- Prevent changing directories when saving files
- Disable delete, rename and right-click in file dialog boxes
- Default configuration
- Directory for spectral data
- Directory for experiments
- Directory for log files
- Directory for quant methods
- Directory for macros
- Directory for user libraries
- Directory for report templates and notebooks
- Directory for report templates and notebooks

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The screenshot displays the Windows Event Viewer interface. The left-hand tree view shows the 'Thermo Electron' folder selected and circled in red. The main pane shows a list of events for 'Thermo Electron' with columns for Type, Date, Time, Source, Category, and Event ID. A red arrow points from the 'Warning' event at 7:16:49 AM (Source: OMNIC) to the 'Event Properties' dialog box. The dialog box shows the following details:

Field	Value
Date	12/21/2004
Time	7:16:49 AM
Type	Warning
User	AMER\mike.gary
Computer	USMAD-L4F58Y51
Source	Omic
Category	Login Audits
Event ID	28

Description:
Logon failure (unknown user name or bad password).
User : mike.gary

Tracks program use and file events, even when OMNIC is not running!

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

TQ Analyst - [New Method]

File Edit View Diagnostics Window Help

Calibrate Quantify Explain Close Performance Index: N/A Previous: N/A Un-calibrated

Description Pathlength Components Standards Regions Corrections Other Report

Suggest... How To...

Revision: 0
 Last saved: Tue Mar 29 15:53:03 2005

Method Title

Analysis Type

Quantitative analysis

- Simple Beer's law
- Classical least squares (CLS)
- Stepwise multiple linear regression (SMLR)
- Partial least squares (PLS)
- Principal component regression (PCR)
- Undecided

Classification

- Similarity match
- Distance match
- Discriminant analysis
- Search standards
- QIC Compare search

Measurement

- Measurement only

Developer's Name

Method Description

Quantitative analysis

- Simple Beer's law
- Classical least squares (CLS)
- Stepwise multiple linear regression (SMLR)
- Partial least squares (PLS)
- Principal component regression (PCR)

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